

Protocol for ER

Obtained in run B31, user 1593

31 Dec 2020

Immunostainer

Type: Leica BOND Max

Primary antibody

Clone: IHC423
Producer: GenomeMe
Product no. / lot no.: IHC423 / CC855527
Diluent: GenomeMe Antibody Diluent
Dilution factor: 1:50
Incubation time / temperature: 30 min. / 20°C

Epitope retrieval, HIER

Device: On Board / On Machine
Buffer: Leica Bond Epitope Retrieval Solution 2
Heating time at max. temp.: 30 min.
Maximum heating temp.: 100°C

Visualization system

Producer: Leica
Product / no: Bond Refine / DS9800
Incubation time linker: 8 min.
Incubation time polymer: 8 min.
Incubation temperature: 20°C

Chromogen

Producer: Leica
Product / no: Bond Refine / DS9800
Incubation time / temperature: 10 min. / 20°C
Enhancement: None